malleonite

REF. 90062

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Night splint

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Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.;

SELECTION/SIZE

Size
Shoe size

S/M	L/XL
36/42	43/48

Color: black
Right and left

CHARACTERISTICS AND MATERIALS

- The hinged polymer splint can be positioned for personalised support using the adjustable plantar flexion strap
- Comfortable, soft padding
- The adjustable Velcro® calf and foot straps are padded for optimum comfort.

INDICATIONS

- Plantar fasciitis
- Achilles tendinopathy
- Foot drop
- Pain upon awakening

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE

English

- 1 Undo all the straps and fasten them back on themselves for the moment.
- 2 Slide the foot into position and rest the splint of the top of the foot and the ankle (fig. A).
- 3 Pull the Velcro® calf strap and fasten (fig. B).
- 4 Next pull the foot (fig. C) strap and fasten.
- 5 [For the first time of wearing only] Adjust the dorsal flexion using the Velcro® strap (fig. D)

MAINTENANCE

do not bleach

do not iron

washing instructions: Hand wash in lukewarm water with mild soap; leave to dry away from direct heat.

no chemical cleaning

do not tumble-dry

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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Headquarter: ORTHOSERVICE AG
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland