linealastik

genulastik 04

malleolastik 37

epilastik 32

manulastik 19

REF.90104 - Genulastik04

REF.90137 - Malleolastik37

REF.90132 - Epilastik32

REF.90119 - Manulastik19

Take
care
feel
better

ORTHO SERVICE
R+T

genulastik 04

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Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

MAINTENANCE

do not bleach

do not iron

washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

no chemical cleaning

do not tumble-dry

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The orthosis should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or (ref.90104) strong electromagnetic fields.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician.

Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

CHARACTERISTICS AND MATERIALS

- Stretch knit fabric featuring areas of different targeted compression
- Anatomically designed
- Silicone patella stabiliser
- Fine side splints keep the knee support firmly in place
- The upper edging is strengthened with silicone trim so the support feels softer and keeps its shape when being pulled on.

INDICATIONS

- Gonalgia from bruising
- Early stage arthrosis
- Femoro-patellar chondropathy and chondromalacia in the knee
- Mild instability of the patella
- Prevention of displacement of the patella during sporting activity.

CONTRAINDICATIONS

Currently no known.

PUTTING ON THE APPLIANCE

- Slide on the knee support with the label facing upwards
- Arrange the knee support so that the patella is correctly positioned within the dedicated silicone stabiliser.

CHARACTERISTICS AND MATERIALS

- Stretch knit fabric featuring areas of different targeted compression
- Anatomically designed
- Silicone inserts shaped to fit the malleoli
- The upper edging is strengthened with silicone trim so the support feels softer and keeps its shape when being pulled on.

INDICATIONS

- Grade 1 ankle sprains
- Tendinopathy of the instep and peroneal tendons.

CONTRAINDICATIONS

Currently no known

PUTTING ON THE APPLIANCE

- Slide on the ankle support with the label facing upwards
- Arrange the ankle support so that the malleoli are correctly positioned within the dedicated silicone inserts.

CHARACTERISTICS AND MATERIALS

- Stretch knit fabric featuring areas of different targeted compression
- Anatomically designed
- Silicone inserts shaped to redistribute compression on the epicondyle and epitrochlea and pressure pads on the forearm
- The edging is strengthened with silicone trim so the support feels softer and keeps its shape when being pulled on.

INDICATIONS

- Lateral epicondylitis (tennis elbow) and medial epicondylitis (golfer's elbow)
- Chronic pain from bursitis or tendinitis.

CONTRAINDICATIONS

Currently no known.

PUTTING ON THE APPLIANCE

- Slide on the elbow support with the label facing upwards: the larger part of the silicone insert should rest on the forearm
- Position the elbow support so that the openings in the silicone insert are centred around the epicondyle and epitrochlea.

CHARACTERISTICS AND MATERIALS

- Stretch knit fabric featuring areas of different targeted compression
- Anatomically designed
- Thumb hole
- Polymeric palm splint
- Fastens with a Velcro® wrist strap
- Silicone insert shaped to fit the styloid process
- The forearm edging is strengthened with silicone trim so the support feels softer and keeps its shape when being pulled on.

INDICATIONS

- Wrist and carpal sprains
- Tendinitis of the wrist
- Arthrosis and arthritis of the wrist
- After-effects of wrist fractures
- Post-surgery wrist treatment.

CONTRAINDICATIONS

Currently no known.

PUTTING ON THE APPLIANCE

- Slide on the wrist support with the label facing upwards
- Position the wrist support so that the thumb goes through the hole provided and the stiff splint rests against the palm of the hand
- Fasten the Velcro® wrist strap.

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Headquarter: ORTHOSERVICE AG
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland