top -S 10°/15°
top II -S 10°/15°
top III -S 10°
top4 -S 10°/15°

REF. E110 10° · REF. E210 15°

REF. E111 10° · REF. E211 15°

REF. E112 10°

REF. E410 10°/15°

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Shoulder support

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Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

REF. E110 10° · REF. E210 15°

Size	S	M	L
Forearm length cm*	< 34	35/39	> 40

*Measurement taken from elbow to end of closed fist

Colour: black
fits right and left

REF. E111 10° · REF. E211 15° · REF. E112 10°

Size	S	L
Forearm length cm*	< 38	>39

*Measurement taken from elbow to end of closed fist

Colour: black
fits right and left

REF. E410 10°/15°

Size	One size
------	----------

Colour: black
Right or left

CHARACTERISTICS AND MATERIALS

- High-density foam covered with a fabric woven with silver fibres, which goes next to the skin.
- The pillow has support at the front for optimum stability; (only REF.410) the armrest area is angled for enhanced comfort.
- A comfortable sleeve protects the neck.
- Brushed nylon on the outside of the cover and sleeves fastening the Velcro®.
- The sling is made from a breathable material.
- Straps with buckles and adjustable Velcro® fasteners.
- Thumb support strap.

INDICATIONS

Post-operative care (in abduction) after operations on the rotator cuff, the joint of the glenoid labrum, capsular stabilisation, operations to the soft tissue of the shoulder

CONTRAINDICATIONS

Currently no known

ADAPTING FOR THE DOCTOR / TECHNICIAN


Top -S 10°/15° · REF. E110 10° · REF. E210 15°

Top II-S 10°/15° · REF. E111 10° · REF. E211 15°

Top 4-S 10°/15° · REF. E410

English

- 1 Undo the two side-release buckles.
(REF. 410 ONLY) If necessary, detach the sling from the pillow and reposition it on the actual pillow so it matches up with the marking for the desired degree of abduction: orange=10°; white=15° (fig. A1).
(EXCEPT FOR REF. 410) If necessary, detach the sling from the pillow and reposition it with Velcro® on the actual pillow on the same side as the shoulder requiring immobilisation. (fig. A2)
- 2 Position the pillow on the hip on the same side as the shoulder requiring immobilisation.
- 3 Fasten the waist strap with the side-release buckle (fig. B) and tighten with the corresponding Velcro® (fig. C). If it is too long, cut to size and reposition the Y-shaped Velcro® strap.
- 4 Carefully slide the forearm into the sling (fig. D1), fasten it with the Velcro strap provided and (REF. 410 ONLY) adjust the length of the back strap (fig. D2).
- 5 Pass the vertical strap over the opposite shoulder and fasten the side-release buckle (fig. E).
- 6 Tighten the shoulder strap with the corresponding Velcro®. If it is too long, cut to size and reposition the Y-shaped Velcro® strap.
- 7 Position the padded sleeve to protect the neck.
- 8 Position the thumb support (and exercise ball) and attach it to the pillow with the Velcro® (fig. F).

PUTTING ON THE APPLIANCE FOR THE PATIENT

Top -S 10°/15° · REF. E110 10° · REF. E210 15°

Top II-S 10°/15° · REF. E111 10° · REF. E211 15°

Top 4-S 10°/15° · REF. E410

English

Maintain the same configuration as the first time the support was worn.

To remove the support open the sling and undo the buckles.

The shoulder support can be redonné as follows:

- 1 Position the pillow on the hip on the same side as the shoulder requiring immobilisation and rest the arm in the sling.
- 2 Fasten the waist strap with the side-release buckle (fig.A).
- 3 Close the sling with the Velcro fasteners and slide the thumb into its support.
- 4 Pass the vertical strap over the opposite shoulder and fasten the side-release buckle (fig.B).

Top III-S 10° · REF. E112

English

- 1 If necessary, remove the arm sling from the cushion and attach to the cushion using the Velcro® fastener so that the orthotic lies to the left or the right-hand side.
- 2 Place the cushion laterally at the side to which the shoulder should be fixed, close the waist belt with the push-in buckle fastener (fig. A) and tighten with the corresponding Velcro® fastener (fig. A). Should the belt be too long, then it can be shortened. To do this, cut the belt to the appropriate length and re-attach the Y-shaped Velcro® fastener end.
- 3 Carefully insert the forearm into the arm fixation (fig. B); then close the arm fixation with the respective Velcro® fasteners at the wrist (fig. C) and the forearm (fig. D).
- 4 Insert the arm on the opposite side from the shoulder to be fixed through the loop created by the shoulder and underarm belts (fig. E) and close the orthotic with the push-in buckle fastener on the arm fixation (fig. F).
- 5 Close all the belts pertaining to the healthy shoulder with the respective Velcro® fasteners (fig. G); Back belt, underarm belt, shoulder belt and front belt. If the belt is too long, it can be shortened; to do this, cut the belt to the appropriate length and re-attach the Y-shaped Velcro® fastener end.
- 6 Check the placement of the underarm and neck padding; re-position the padding if required (fig. H).
- 7 Attach the thumb support (and the softball for rehabilitation) onto the cushion with the Velcro® fastener (fig. I).

Top III-S 10° · REF. E112

Retain the configuration as adjusted by the doctor/technician. In order to remove the orthotic, open the arm fixation and unhook the push-in buckle fasteners. To put the orthotic on again, proceed as follows:

- 1 Place the cushion laterally at the side on which the shoulder is to be fixed; then place the arm into the breathable arm fixation.
- 2 Close the waist belt with the push-in buckle fastener (fig. A).
- 3 Close the arm fixation with the respective Velcro® fasteners and insert the thumb into the thumb support.
- 4 Insert the arm on the opposite side from the shoulder to be fixed through the loop created by the shoulder and underarm belts (fig. B); Close the orthotic with the push-in buckle fastener on the arm fixation (fig. C).

MAINTENANCE

do not bleach

do not iron

washing instructions: Wash the cushion cover (after the foam inner padding has been removed), belts and the neck and underarm padding by hand with warm water (at max. 30°C) and neutral soap; rinse thoroughly. Dry away from heat sources.

no chemical cleaning

do not tumble-dry

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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**Headquarter: ORTHOSERVICE AG**
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Orthoservice AG (CH) · 6830 Chiasso (TI) · Switzerland